

Experiment 4: The Qual Scheme


Selected Anion Analysis

Anion Elimination Tests



In this experiment you will analyze for the presence of five anions. The knowns are provided as sodium salts (all white crystals) and 0.2M solutions of the same salts. The anion unknown provided is a mixture of the solids ground together. While many anions are known by their intrinsic color, none of these five are colored. Nor will there be any keys to identification by knowing the cation since all sodium salts are soluble. Known samples are to be used separately in three experiments described below. In these experiments you may want to join with a partner.

Experiment 1 (Essentially Procedure 22)



Using a few drops of each of the 0.2M solutions of anion knowns in separate labeled tubes, add 0.2M AgNO_3 . Note which anions produce precipitates (all are insoluble silver salts) and note the color of each. For each anion write an equation of this form:



Anion Experiment One Diagram




Anion Experiment Two Diagram




Experiment 2 (Essentially Procedure 23)

Using two drops of each of the known solutions in separate labeled tubes, first add a drop of 5M NH₃ to each. Then add 2 or 3 drops of 0.2M BaCl₂. Note which produce precipitates and their color. Write equations for precipitate formations in this form:



Where precipitates are formed, centrifuge and discard the decantates. Treat the precipitates with 2M HCl (look for bubbles). Write equations for precipitate dissolution in acid in this form:



Fill out the flow diagram “Anion Experiment Two” with the formula of products and colors on the basis of your experimental results.

Experiment 3 (Essentially Procedure 24)



Take a small pinch of each of the five known sodium salts of the anions in five separate dry test tubes and treat each successively with 2 or 3 drops of 18M H₂SO₄. At this time let us remind you how corrosive is concentrated H₂SO₄. Treat it with great respect. Don't point the tube towards anyone, especially later when heating. And of course many acid gases are produced; smell them only by wafting. As you add the conc H₂SO₄, watch for effervescence (bubbles) and discoloration in the acid solution.

Two of the five salts produce gases. Carbonates produce the mildly acid, colorless, odorless gas CO₂, with any strong acid solution. But the other will bubble only with very concentrated strong acids. Chlorides, on treatment with 18 M H₂SO₄ produce HCl, a colorless sharply acidic gas which fumes in moist air.

Write two equations for the formation of these gases.



Then, with the three salts which did not produce gases, take each mixture and heat it carefully till hot and look for further activity. You will note that nitrate produces yet another yellow brown gas by decomposition of the nitric acid formed.



Procedure 22

When analyzing the “anion unknown,” remember that you must interpret the results in view of the fact that more than one anion is present. In the case of ‘single salt’ unknown consider whether any acid was used to put the salt in solution.

Place three drops of the solution to be analyzed in a test tube and add 2 drops of 0.2M AgNO₃. Note the color of any precipitate. Centrifuge and test for complete precipitation with 0.2M AgNO₃. Centrifuge and discard the decantate. Treat the precipitate with plenty of 3M HNO₃, look for bubbles and any remaining precipitate. Centrifuge and treat the decantate with 15M NH₃ dropwise until just basic, then with 5M HC₂H₃O₂ until just acid. Add another drop of 0.2M AgNO₃ and note the color of any precipitate. Under these conditions no problem with AgOH or Ag₂O should occur.

Procedure 23

As with Proc 22 interpretation will depend on whether the sample has several or a single anion. To three drops of solution to be analyzed add 5M NH₃ until the sample is just basic, usually one drop if the sample had not been acidified. Add 2 drops of 0.2M BaCl₂. If a precipitate is produced, isolate the precipitate and treat it with 6M HCl. Look for bubbles and note if any precipitate remains.

Procedure 24

Treat a small pinch of the solid unknown to be tested with 3 drops of 18M H₂SO₄. Look for bubbles, color in the solution, color in the gas, and test the gas with wet blue litmus. If strongly acid, colored gases are not observed, carefully heat the mixture over a small flame. Now yellow-brown gas on heating indicates nitrate. These latter tests are usually feasible only if no strongly acidic gas producing anions are present.

Specific Anion Procedures

Procedure 25 Test for Borate Ion

With a small quantity of the solid salt to be tested for borate in a casserole, add 3 drops 18M H₂SO₄ (hood) and 20 drops of methanol. Heat this mixture over a burner in the hood about 5 seconds then set it afire by tipping the casserole so the flame lights the fumes. The immediate appearance of a green flame indicates the presence of borate ion. Though the cations Ba²⁺ and Cu²⁺ have distinct green flames they are not volatile and don't give the immediate color.



Equations for Procedure 25



Procedure 26 Test for Carbonate Ion

Place a small amount of solid sample of the salt or mixture to be tested in a dry test tube. Get about one mL of saturated Ba(OH)₂ from the reagent shelf. If this solution is not clear, centrifuge it and draw some up in your Pasteur pipet. Now, squirt a portion of 2M HCl onto the solid you are testing. If carbonate is present there will usually be vigorous effervescence. Insert the Pasteur pipet about half way into the tube and squeeze out a drop so that it hangs on the pipet. Immediate clouding of the suspended drop implies the presence of CO₂, thus carbonate. If any Ba(OH)₂ is dropped in the solution it doesn't matter.

Equations for Procedure 26



white

Procedure 27 Test for Sulfate Ion

Use a few drops of the solution to be tested and acidify it with 6M HCl. Add one or two drops of 0.2M BaCl₂. A white precipitate indicates the presence of sulfate ion.

Equation for Procedure 27